


### Post Mortem Needles 3/8 Circle Cutting Edge

**08816:** non-sterile, stainless  
12/pkg

Available Sizes	Dimensions
3	1.42mm x 75mm
4	1.80mm x 100mm
5	2.00mm x 125mm
6	2.00mm x 152mm



### Straight Cutting Needles Triangular Point

**08817:** non-sterile, stainless  
12/pkg


Available Sizes	Dimensions
2	1.02mm x 2"
3	1.27mm x 3"
4	1.42mm x 4"



### Post Mortem Needles Serpentine Double Curved Cutting Edge

**08819:** non-sterile, stainless  
12/pkg


Available Sizes	Dimensions
0	5"
1	4.75"
2	4.25"
3	3.25"
4	2.25"



### Post Mortem Needles Half Curved Cutting Edge

**08820:** non-sterile, stainless  
12/pkg

Available Sizes	Dimensions
00	2.20mm x 150mm
0	2.20mm x 140mm
1	2.00mm x 114mm
2	1.8mm x 89mm




All non-sterile  
needles are sold  
in multiples of  
10 dozen.

\*Needles are not shown to actual size.

## Regular Surgeon Needles Half Curved Cutting Edge

**08821:** non-sterile, stainless  
12/pkg


Available Sizes	Dimensions
1	1.21mm x 80mm
2	1.20mm x 68mm
4	1.20mm x 64mm
6	1.10mm x 55mm
8	1.00mm x 50mm
10	.93mm x 43mm
12	.76mm x 35mm
14	.69mm x 29mm
16	.60mm x 27mm
18	.64mm x 25mm
20	.53mm x 22mm
22	.50mm x 19mm



## Regular Surgeon Needles 3/8 Circle Cutting Edge

**08822:** non-sterile, stainless  
12/pkg


Available Sizes	Dimensions
1	1.20mm x 80mm
2	1.22mm x 68mm
4	1.21mm x 64mm
6	1.10mm x 55mm
8	1.00mm x 50mm
10	.94mm x 43mm
12	.76mm x 35mm
14	.69mm x 29mm
16	.64mm x 27mm
18	.53mm x 22mm
20	.48mm x 17mm



## Regular Surgeon Needles 3/8 Circle Taper Point

**08822T:** non-sterile, stainless  
12/pkg

Available Sizes	Dimensions
4	1.16mm x 65mm
6	1.09mm x 55mm
12	.76mm x 35mm
14	.69mm x 29mm
16	.64mm x 27mm
18	.53mm x 22mm




## Mayo's Catgut Needles 1/2 Circle Taper Point

**08824:** non-sterile, stainless  
12/pkg

Available Sizes	Dimensions
2	1.40mm x 45mm
3	1.40mm x 38mm
4	1.40mm x 35mm
6	1.30mm x 28mm

**09824:** sterile, stainless  
50/box


Available Sizes	Dimensions
2	1.40mm x 45mm
3	1.40mm x 38mm
4	1.40mm x 35mm
5	1.30mm x 32mm
6	1.30mm x 28mm



## Mayo's Catgut Needles 1/2 Circle Trocar Point

**09826:** sterile, stainless  
50/box

Available Sizes	Dimensions
2	1.40mm x 45mm
3	1.40mm x 38mm
4	1.40mm x 35mm
5	1.30mm x 32mm
6	1.30mm x 28mm




## Keith's Straight Abdominal Needles Triangular Point

**08827:** non-sterile, stainless  
12/pkg

Available Sizes	Dimensions
1 1/2"	.74mm x 38mm
1 3/4"	.75mm x 44mm
2	.76mm x 51mm
2 1/2"	.81mm x 64mm
2 7/8"	.92mm x 76mm
3 7/8"	1.00mm x 95mm

**09827:** sterile, stainless  
50/box

Available Sizes	Dimensions
2	.76mm x 51mm
2 1/2"	.81mm x 64mm
2 7/8"	.83mm x 76mm



All non-sterile  
needles are sold  
in multiples of  
10 dozen.

\*Needles are not shown to actual size.

**Eye Needles**

**3/8 Circle Cutting Edge**

**08830:** non-sterile, stainless

12/pkg

Available Sizes	Dimensions
1	.53mm x 28mm
2	.53mm x 22mm



All non-sterile needles are sold in multiples of 10 dozen.

**Fistula Needles**

**1/2 Circle Cutting Edge**

**08832:** non-sterile, stainless

12/pkg

Available Sizes	Dimensions
0	1.42mm x 64mm
9	.76mm x 19mm




**Regular Surgeon Needles**

**1/2 Circle Cutting Edge**

**08834:** non-sterile, stainless

12/pkg

Available Sizes	Dimensions
00	1.20mm x 68mm
0	1.20mm x 64mm
1	1.09mm x 59mm
2	1.00mm x 50mm
4	.90mm x 43mm
6	.68mm x 29mm
7	.60mm x 27mm
8	.55mm x 21mm
9	.50mm x 19mm




**Regular Surgeon Needles**

**1/2 Circle Taper Point**

**08834T:** non-sterile, stainless

12/pkg

Available Sizes	Dimensions
0	1.20mm x 64mm
2	1.00mm x 50mm
4	.90mm x 43mm
6	.68mm x 29mm
9	.53mm x 19mm



**Mayo's Intestinal Needles**

**3/8 Circle Taper Point**

**08837:** non-sterile, stainless

12/pkg

Available Sizes	Dimensions
1	.53mm x 25mm




**Ferguson's Needles**

**1/2 Circle Taper Point**

**08842:** non-sterile, stainless

12/pkg

Available Sizes	Dimensions
2	1.27mm x 70mm
4	1.22mm x 60mm
6	1.10mm x 55mm
8	.92mm x 44mm
10	.92mm x 39mm
12	.76mm x 35mm
14	.76mm x 29mm
16	.70mm x 25mm
18	.53mm x 19mm
20	.64mm x 15mm



**Lane's Cleft Palate Needles**

**1/2 Circle Cutting Edge**

**08844:** non-sterile, stainless

12/pkg

Available Sizes	Dimensions
1	.71mm x 25mm




\*Needles are not shown to actual size.

## Loopuyt's Cutting Edge Needles

**08846:** non-sterile, stainless

12/pkg

Available Sizes	Dimensions
1	1.42mm x 80mm
2	1.22mm x 68mm
3	1.10mm x 60mm
4	1.00mm x 52mm
5	1.00mm x 46mm



All non-sterile needles are sold in multiples of 10 dozen.

## Lane's Cleft Palate Needles

### 3/8 Circle Cutting Edge

**08848:** non-sterile, stainless

12/pkg

Available Sizes	Dimensions
1	.71mm x 22mm
2	.69mm x 19mm
2	.69mm x 16mm



## Davis Tonsil Needles

### 1/2 Circle Taper Point

**09849:** sterile, stainless

50/box

Available Sizes	Dimensions
one size	1.20mm x 25mm



## Martin's Uterine Needles

### 1/2 Circle Cutting Edge

**08860:** non-sterile, stainless

12/pkg

Available Sizes	Dimensions
1	1.60mm x 74mm




## Regular Surgeon Needles

### 3/8 Circle Reverse Cutting Edge

**08922:** non-sterile, stainless

12/pkg

Available Sizes	Dimensions
6	1.10mm x 55mm
8	1.00mm x 50mm
10	.93mm x 43mm
12	.76mm x 35mm
14	.69mm x 29mm
16	.64mm x 27mm
18	.53mm x 22mm



## Regular Surgeon Needles

### 1/2 Circle Reverse Cutting Edge

**08934:** non-sterile, stainless

12/pkg

Available Sizes	Dimensions
3	1.00mm x 50mm
5	.85mm x 38mm
7	.60mm x 27mm

